## **Protein Solutions™**

## utomated, Online Second Virial Coefficient (A2) Measuremen

Provide the second virial coefficient,  $(A_2)$ .

 $A_2$  values are dependent on both the protein and its solvent: Temperature, salinity, pH, chemical excipients, etc., can induce dramatic shifts in  $A_2$ . Positive  $A_2$ corresponds to repulsive intermolecular interaction, whereas negative  $A_2$  corresponds to attractive intermolecular interaction. Rapid determination of the buffer conditions yielding an optimal  $A_2$ , such as a positive  $A_2$ for formulations, or a slightly negative  $A_2$  for crystallization, is highly desirable.

In the past, such measurements were labor intensive, time consuming, and expensive in terms of time and protein: For each buffer condition tested, extensive dialysis and preparation of multiple concentrations were required. Additionally, each concentration had to be injected manually. Wyatt Technology has overcome these issues, and developed an accurate, automated method for measuring  $A_2$  online. In our online  $A_2$  method, only a single protein stock solution is needed for an entire series of measurements.

In addition to complete automation and online dialysis, the online  $A_2$  method has the added benefit of flushing the light scattering flow cell between each injection, helpful for researchers working with "sticky" proteins.

The Trainoff-Wyatt Online  $A_2$  method is as follows: A sample delivery unit, such as a standard autosampler, delivers a sequence of 5 or more different volumes (such as 20 $\mu$  to 200 $\mu$ L) through an optional size exclusion chromatography column (SEC) or a desalting column, to a DAWN or miniDAWN light scattering instrument and an Optilab RI. The Optilab RI serves as a concentration detector and also records the complete dialysis of each injection.

The Online  $A_2$  method will enable researchers to fully utilize the second viral coefficient as an important indicator of protein-protein interaction. The historic impediments to this type of measurement, lengthy sample preparation and experimental labor, have been overcome. Now, using standard chromatography equipment and Wyatt detectors, an entire series of  $A_2$ values can be determined in less then a day!



Figure 1: Upper right graph (foreground) shows the light scattering signal (Detector 11) of the  $20\mu$ L to  $200\mu$ L antibody solution injection series. The expanded (box) region shows both the light scattering signal in black, and differential refractive index signal in red, demonstrating complete dialysis of injected antibody with the mobile phase by the SEC column. Only the protein peaks are used in the Online A, analysis.



Figure 1: Zimm Plot from the light scattering and concentration data of the protein peaks. Each color represents a different peak. Note the excellent agreement of the measured data to the fitted data. This Zimm plot reports the antibody in PBS has an  $A_2$ of -5.46E-5 mol mL/g2.

