




Who knew you had the budget for ten extra people.



Now you can add staff without staffing up. The DynaPro™ Plate Reader makes you at least 10 times more productive than you are with a batch dynamic light scattering instrument. On a single well plate – with as few as 5 microliters of sample – you can measure the sizes of up to 1,536 proteins, liposomes, biomacromolecules, viral or nanoparticles without ever having to change or clean a cuvette again. It's also temperature controlled to help you perform a variety of experiments faster than ever. In fact, the DynaPro Plate Reader is so productive it pays for itself in just a few months. And it interfaces with liquid handling robots to speed things up even more. For more information on how to run more samples without running around, call **805.681.9009** or visit www.wyatt.com

