



**WYATT**  
TECHNOLOGY

## **NaCl Solution Standards**

6330 Hollister Avenue  
Santa Barbara, CA 93117

Tel: +1 (805) 681-9009  
Fax: +1 (805) 681-0123  
Web: www.wyatt.com

# Certificate of Analysis

## Sodium Chloride Solutions Kit

**Wyatt p/n: P8400-B**  
**Lot #: 02062019**

**Catalog No.: P8400**  
Storage: Room Temperature  
Expires: January 31, 2022

### Standard Analysis Data

Standard	CAS#	Percent Purity	Concentration (weight/volume)
18 mega-ohm Water	7732-18-5	100.00%	NEAT
0.100 mg/mL NaCl Solution	7647-18-5	99.999%	0.100 mg/mL
0.250 mg/mL NaCl Solution	7647-18-5	99.999%	0.250 mg/mL
0.500 mg/mL NaCl Solution	7647-18-5	99.999%	0.500 mg/mL
1.000 mg/mL NaCl Solution	7647-18-5	99.999%	1.000 mg/mL
1.200 mg/mL NaCl Solution	7647-18-5	99.999%	1.200 mg/mL
2.000 mg/mL NaCl Solution	7647-18-5	99.999%	2.000 mg/mL
3.000 mg/mL NaCl Solution	7647-18-5	99.999%	3.000 mg/mL
4.000 mg/mL NaCl Solution	7647-18-5	99.999%	4.000 mg/mL
5.000 mg/mL NaCl Solution	7647-18-5	99.999%	5.000 mg/mL

**$dn/dc$  at 658 nm: 0.174 mL/g**  
 **$dn/dc$  at 690 nm: 0.172 mL/g**

Certification: Michelle Chen

Date: March 26, 2019

Dr. Michelle Chen, Director of Analytical Services

[P1154-B