



WYATT
TECHNOLOGY

6300 Hollister Avenue
Santa Barbara, CA 93117

Tel: +1 (805) 681-9009
Fax: +1 (805) 681-0123
Web: www.wyatt.com

NaCl Solution Standards

Certificate of Analysis

Sodium Chloride Solutions Kit

Wyatt p/n: P8400

Lot #: 082812

Catalog No.: P8400

Storage: Room Temperature

Expires: September 2017

Standard Analysis Data

Standards	CAS#	Percent Purity	Concentration (weight/volume)
18 mega-ohm Water	7732-18-5	100.00%	NEAT
0.100 mg/mL NaCl Solution	7647-18-5	99.999%	0.100 mg/mL
0.250 mg/mL NaCl Solution	7647-18-5	99.999%	0.250 mg/mL
0.500 mg/mL NaCl Solution	7647-18-5	99.999%	0.500 mg/mL
1.000 mg/mL NaCl Solution	7647-18-5	99.999%	1.000 mg/mL
1.200 mg/mL NaCl Solution	7647-18-5	99.999%	1.200 mg/mL
2.000 mg/mL NaCl Solution	7647-18-5	99.999%	2.000 mg/mL
3.000 mg/mL NaCl Solution	7647-18-5	99.999%	3.000 mg/mL
4.000 mg/mL NaCl Solution	7647-18-5	99.999%	4.000 mg/mL
5.000 mg/mL NaCl Solution	7647-18-5	99.999%	5.000 mg/mL

dn/dc at 658 nm: 0.174 mL/g

dn/dc at 690 nm: 0.172 mL/g

Certification: ***Michelle Chen***

Dr. Michelle Chen, Director of Analytical Services

Date: September 5, 2012