

**Wyatt
Technology**
CORPORATION

6300 Hollister Avenue • Santa Barbara, CA 93117
TEL (805) 681-9009 • FAX (805) 681-0123
E-mail: info@wyatt.com • URL <http://www.wyatt.com>

NaCl Standards

Certificate of Analysis Sodium Chloride Solutions Kit


Wyatt p/n: P8400
External Lot #: 102510

Catalog No. P8400
Storage: Room Temperature

Standard Characterization Data

Components	Lot #	CAS#	Percent Purity	Concentration (weight/volume)
18 mega OHM Water	A077357	7732-18-5	100.00%	NEAT
0.1000 mg/ml NaCl Solution	A077110	7647-14-5	99.999%	1.0000e-04 g/ml
0.2500 mg/ml NaCl Solution	A077114	7647-14-5	99.999%	2.5000e-04 g/ml
0.5000 mg/ml NaCl Solution	A077118	7647-14-5	99.999%	5.0000e-04 g/ml
1.0000 mg/ml NaCl Solution	A077195	7647-14-5	99.999%	1.0000e-03 g/ml
1.2000 mg/ml NaCl Solution	A077199	7647-14-5	99.999%	1.2000e-03 g/ml
2.0000 mg/ml NaCl Solution	A077205	7647-14-5	99.999%	2.0000e-03 g/ml
3.0000 mg/ml NaCl Solution	A077249	7647-14-5	99.999%	3.0000e-03 g/ml
4.0000 mg/ml NaCl Solution	A077251	7647-14-5	99.999%	4.0000e-03 g/ml
5.0000 mg/ml NaCl Solution	A077256	7647-14-5	99.999%	5.0000e-03 g/ml

dn/dc at 658 nm: 0.174 ml/g
dn/dc at 690 nm: 0.172 ml/g



**Wyatt
Technology**
CORPORATION

6300 Hollister Avenue • Santa Barbara, CA 93117
TEL (805) 681-9009 • FAX (805) 681-0123
E-mail: info@wyatt.com • URL <http://www.wyatt.com>

NaCl Standards

Certificate of Analysis Sodium Chloride Solutions Kit

Wyatt p/n: P8400
External Lot #: 102510

Catalog No. P8400
Storage: Room Temperature

Standard Characterization Data

Components	Lot #	CAS#	Percent Purity	Concentration (weight/volume)
18 mega OHM Water	A077357	7732-18-5	100.00%	NEAT
0.1000 mg/ml NaCl Solution	A077109	7647-14-5	99.999%	1.0000e-04 g/ml
0.2500 mg/ml NaCl Solution	A077113	7647-14-5	99.999%	2.5000e-04 g/ml
0.5000 mg/ml NaCl Solution	A077117	7647-14-5	99.999%	5.0000e-04 g/ml
1.0000 mg/ml NaCl Solution	A077195	7647-14-5	99.999%	1.0000e-03 g/ml
1.2000 mg/ml NaCl Solution	A077199	7647-14-5	99.999%	1.2000e-03 g/ml
2.0000 mg/ml NaCl Solution	A077202	7647-14-5	99.999%	2.0000e-03 g/ml
3.0000 mg/ml NaCl Solution	A077249	7647-14-5	99.999%	3.0000e-03 g/ml
4.0000 mg/ml NaCl Solution	A077253	7647-14-5	99.999%	4.0000e-03 g/ml
5.0000 mg/ml NaCl Solution	A077256	7647-14-5	99.999%	5.0000e-03 g/ml

dn/dc at 658 nm: 0.174 ml/g
dn/dc at 690 nm: 0.172 ml/g